



PRECISION LAB FURNITURE INDUSTRIES



MODULAR WORK TABLES

Island table forms the central working area of any laboratory. The table made of high grade galvanized steel finished in epoxy polyester powder coating and high quality black granite working top. Island table with top drawer and bottom cupboard of removable and accessible type on both sides with knee space ensure maximum working and storage space





FEATURES

One Reagent rack of standard size made out of 2 tier steel partition with tray. Two 15/5 amps 3 pin electrical socket. Polypropylene sink with 3 way water tap

WALL STORAGE CABINET

Wall Table consists of individual modular with top drawer bottom cup board arrangement. The cupboard is provided with adjustable horizontal partition. The working surface will be granite of uniform thickness. The table will be provided with electrical panel upon customer request.



LABORATORY FUME EXHAUST HOODS

Laboratory Fume Exhaust hoods are available in metal construction with inner FRP inner or in complete FRP constructions. Heavy-duty and dynamically balanced blower assembly provided in the unit removes smoke, fume, odor from any chemical process conducted with finely polished granite. Provided suitable lighting inside the work chamber. Available with or without storage cupboard under the workbench available in ductless type.



CHEMICAL STORAGE CUPBOARD

Chemical storage cabinet provided to store solvent & chemical inside laboratory. No. of racks are provided worth front of glass door. Exhaust port on top allow no. of cabinet to be connected to common blower in cluster, size can be altered as per requirement.



TISSUE CULTURE RACK

The racks consists of Epoxy powder coated M.S Tubular frame work fitted with horizontal racks (5 to 6 Nos required). Two to four sets of fluorescent tube light will be provided to give the required illumination lux. A control panel and timer will be provided on the side in incorporating the individual chokes and switches for the illumination. The rack will have castor wheel for easy mobility.



S.S LAB STOOL.

DIMENSSION	: 300dia x 600mm H (750mm)
TYPE	: Revolving height adjustment
M.O.C	: Stainless steel / GI Powder Coated

ANTI -VIBRATION TABLE

Anti - vibration table used with highly sensitivity instruments, like electro balance, galvanometer, and calibration work etc. the outer table is positioned independently of the inner table. The inner table is a heavy weighted pre-fabricated. System which is mounted on specially rubber mounts to dampen the minutes vibration.



EYE WASHER


Eye wash is a unit that is used to irrigate and flush the eyes and cover the entire body.



Table Top



Island Top



L Top



U Top



Linear Top

Modular Cabinet




1 Draw 2 Door




1 Draw 1 Door




3 Draw



2 Door



1 Door




Sliding Door


Computer Table



Single




Linear



Double Side


Chemical Storage Cupboard




Chemical Storage - 2 Door



Chemical Storage 3 Door



wall Storage - 2 Door




Wall Storage - Sliding Door


Reagent Rack



2 Tier



3 Tier



Continuous


Accessories



3 Way Tap



Single Way Tap




PP Sink



Power Panel for Wall Table



Table Side Panel



Double Side Power Panel



PRECISION LAB FURNITURE INDUSTRIES

503, Veera Boyar Colony, Dr. Nanjappa Road, Coimbatore - 18, TAMILNADU

0422 - 2235970, 2237280, 4506280, 97860 22732, 96269 22732

